
# EXAJET®

- EN** HYDROPHILIC VINYL POLYSILOXANE IMPRESSION MATERIAL IN FOIL POUCH  
**DE** HYDROPHILES VINYL-POLYSILOXAN ABFORMMATERIAL IN SCHLAUCHBEUTELN  
**FR** MATÉRIAU HYDROPHILIQUE POUR EMPREINTES À BASE DE VINYL POLYSILOXANE EN RECHARGE  
**IT** MATERIALE PER IMPRONTE IDROFILICO A BASE DI VINIL POLISILOSSANO IN CONFEZIONE PER DISPENSERS AUTOMATICI  
**ES** MATERIAL DE IMPRESIÓN HIDROFÍLICO DE VINIL POLISILOXANO EN BOLSITA METALIZADA  
**NL** HYDROFIEL VINYLPOLISILOXAAN AFDRUKMATERIAAL IN FOLIEVERPAKKING  
**DA** HYDROFILT A-SILIKONE AFTRYKSMATERIALE I FOLIEPAKKE  
**SV** A-SILIKONBASERAT HYDROFILT AVTRYCKSMATERIAL I FOLIEFÖRPACKNING  
**PT** MATERIAL DE IMPRESSÃO DE VINIL POLISSILOXANO HIDRÓFILO EM BOLSA DE ALUMÍNIO  
**EL** ΥΔΡΟΦΙΛΙΚΗ ΒΙΝΥΛ ΠΟΛΥΣΙΛΟΞΑΝΗ ΑΠΟΤΥΠΩΣΗΣ ΣΕ ΣΥΣΚΕΥΑΣΙΑ ΑΛΟΥΜΙΝΙΟΥ

The logo for GC, consisting of the letters 'G', 'C', and 'I' in a stylized, italicized font. The 'G' and 'C' are connected, and the 'I' is separate. There are small slanted lines above and below the letters.

EN

Prior to use, carefully read the instructions.  
For use only by a dental professional in the recommended indications.

## EXAJET

Hydrophilic vinyl polysiloxane impression material in foil pouch

### RECOMMENDED INDICATIONS

Tray material for both double mix and putty wash impressions techniques. To be used in a dynamic mixing device.

### CONTRAINDICATIONS

In susceptible individuals, sensitisation to the product cannot be excluded. Use of the product should be discontinued if allergic reactions are observed.

### PROPERTIES

	EXAJET Heavy Body Type 1: High consistency Normal set	EXAJET Heavy Body Type 1: High consistency Fast set	EXAMIX Regular Type 2: Medium consistency	EXAMIX Injection Type 3: Low consistency
Total working time (min.sec.)	2'00"	1'30"	2'00"	2'15"
Setting time (min. sec.)	4'00"	3'15"	4'00"	4'00"
Minimum time in mouth (min. sec.)	3'00"	2'30"	4'00"	4'00"
Recovery from deformation (%)	99.4	99.4	99.7	99.7
Maximum strain in compression (%)	2.5	2.5	6.0	7.0
Linear dimensional change after 24 hrs. (%)	0.2	0.2	0.2	0.2

Test conditions : Temperature (23 +/- 2°C). Relative humidity (50 +/- 5%) (ISO 4823 : 1992(E))

### DIRECTIONS FOR USE

#### CARTRIDGE AND MIXING TIP LOADING

#### 1. FILLING THE CARTRIDGE

Insert the foil pouches (1 base and 1 catalyst) with pre-mounted activation caps into the respective chambers of the cartridge. Make sure to align the notches of the cartridge and the activation cap.

## **2. ACTIVATION**

Press the activation caps firmly into the respective chambers of the cartridge until they lock into position.

## **3. USAGE**

1. Insert filled cartridge into a dynamic mixing device.
2. Before attaching a mixing tip, extrude a small amount of material to ensure that it is flowing uniformly from both openings.
3. Replace a used mixing tip with a new one prior to next use.

## **IMPRESSION MAKING**

### **PUTTY WASH TECHNIQUE**

1. Coat the impression tray (e.g. Coe trays) with a thin layer of adhesive and allow to dry.
2. Load EXAJET material into the tray.
3. Take the pre-impresion.
4. Remove after the material has set, clean thoroughly and dry.
5. Prepare the pre-impresion for final impresion.
6. Inject EXAMIX Injection or Regular onto the prepared teeth.
7. Seat tray with pre-impresion into the mouth.
8. Once the setting time has elapsed,remove the impresion from the mouth
9. The resulting impresion should be cleaned and then disinfected according to the manufacturer's recommendations.

### **DOUBLE MIX TECHNIQUE**

1. Coat the impression tray (e.g. Coe trays) with a thin layer of adhesive and allow to dry.
2. Inject EXAMIX Injection or Regular onto the prepared teeth.
3. Load EXAJET material into the tray.
4. Seat tray into the mouth.
5. Once the setting time has elapsed,remove the impresion from the mouth
6. The resulting impresion should be cleaned and then disinfected according to the manufacturer's recommendations.

Note : Make sure to complete step 2 when the tray is fully loaded.

### **NOTES :**

Wait one hour before pouring the model.

EXAJET impresions can be silver- or copper-plated.

**CAUTION :**

1. When extruding EXAJET, care should be taken to avoid mixing or contact with the following materials. They may delay the setting time.
  - Condensation cured silicone impression materials
  - Polysulfide impression materials
  - Polyether impression materials
  - Eugenol containing materials
  - Certain types of latex gloves or retraction cords
  - Also avoid moisture and glycerol when extruding.
2. Care should be taken to avoid getting material on clothing. It is hard to remove when it sets up on clothes.

**STORAGE :**

Store at normal temperature 15-25°C (59-77°F).  
(2-year shelf life guarantee).

**PACKAGES :**

1. Intro package Heavy Body - Normal or Fast set
  - 1 foil pouch base (300 ml) with activation cap
  - 1 foil pouch catalyst (62 ml) with activation cap
  - 1 cartridge
  - 10 dynamic mixing tips
2. Refill package Heavy Body - Normal or Fast set
  - 2 foil pouch base (300 ml each) with activation cap
  - 2 foil pouch catalyst (62 ml each) with activation cap
3. Clinic package Heavy Body - Normal or Fast set
  - 6 foil pouch base (300 ml each) with activation cap
  - 6 foil pouch catalyst (62 ml each) with activation cap
4. Cartridge
5. 50 dynamic mixing tips

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3737 West 127th Street, Alsip, ILL. 60803 U.S.A.
GC ASIA DENTAL PTE. LTD.