

Let's Join Forces!



Smart solutions
from GC
for all your
luting challenges

The logo for GC, featuring the letters 'GC' in a bold, teal, italicized font, with a teal dot above the 'G' and another teal dot above the 'C'.

How to choose the optimal luting cement?

Today more than ever you are spoiled for choice when planning and delivering indirect restoration. From the moment your patient comes in, you have much to consider before deciding on the treatment plan.

As a dentist you know that no two patients are alike and it is therefore important to follow a holistic approach in treatment planning considering the choice of indirect restoration and type of preparation which is best suited to the patient's clinical situation.

An important part of this process is choosing the right luting cement



The key questions you should ask before making the final choice are:
Which is the optimal luting cement to help me

Protect the
tooth
structure?

Protect the
prosthetic
restoration?

The selection of the optimal luting cement is the guarantee for a beautiful smile and long-term patient satisfaction



Discover our High 5 solutions for all your luting challenges

With more than 95 years of experience and one of the most comprehensive product portfolios, we at GC believe we are the ideal partner to guide you towards the most appropriate luting solution.

Because after all, efficient luting is a team effort.

Each of the GC luting solution is suited to different clinical situations and offers you their unique features like physical properties, working time, protocol of use, fluoride release and cost...

Let's Join Forces!

G-CEM LinkAce •
G-CEM Capsule •



• G-CEM LinkForce (dual-cure)
• G-ænial Universal Flo (light-cure)

**HIGH
5
GC LUTING
SOLUTIONS**



• Fuji PLUS
• FujiCEM 2



Freegenol •



• Fuji I
• Fuji TEMP LT



Depending on the clinical situation, the first question you should ask yourself is: Should I CEMENT or BOND?

When to cement?

Cementation is typically used when your preparations allow mechanical retention.
The following clinical aspects will define your choice for the optimal CEMENTATION procedure:

Retentive Preparation



Select **GI** & **RMGI** cement that offers the best balance of ease of use, speed and tooth protection.

Isolation not possible



Select **GI** or **RMGI** cement as they offer moisture tolerance with excellent chemical adhesion to the tooth.

Moderate aesthetics needed



Select **GI** or **RMGI** cement when your clinical situation and/or selected substrate (zirconia, metal & porcelain fused to metal) does not require a very high aesthetic cement.

Patients with special needs



Select **GI** or **RMGI** cement for high caries risk and geriatric patients.

Your options

Conventional glass ionomer (GI) cement | GC's solution: **Fuji I**
Resin-modified glass ionomer (RMGI) cement | GC's solutions: **Fuji PLUS** and **FujiCEM 2**

Source: Dr. Lucile Dahan, France (1) and Prof Jean-François Lasserre, France (2, 3 and 4)

"Most of the time the margins of the preparation are subgingival and total control of humidity proves to be difficult!"

DR. LUCILE DAHAN, FRANCE

When to bond?

Bonding is typically used when extra adhesion is required or when aesthetics are of the utmost importance.

Preparation is not retentive



Select **adhesive resin cement**, when the preparation is not retentive.

Self-adhesive cement (ease of use) can be used for moderate retentive preparations.

Isolation possible



Select **self-adhesive** or **adhesive resin cement** when a dry working field can be assured.

High aesthetic needed



Select **self-adhesive** or **adhesive resin cement** for highly aesthetic outcomes.

Extra adhesion needed



Select **adhesive resin cement** for bonding veneers and hybrid ceramics (e.g. CERASMART).

Your options

Self-adhesive resin cement | GC's solutions: **G-CEM LinkAce** and **G-CEM Capsule**

Adhesive resin cement | GC's solutions: **G-CEM LinkForce (dual-cure)** and **G-ænial Universal Flo (light-cure)**

Discover all of our luting solutions!

With one of the most comprehensive luting portfolios, we at GC can help you maintain your patient's smile.



Our solutions for permanent cementation

Glass ionomer based cements		Composite resin based cements	
<p>Fuji I </p> <p>Conventional glass ionomer (GI) cement</p> <ul style="list-style-type: none"> • Cost-effective solution • Resin-free cement <p>Page 7</p>		<p>G-CEM LinkAce & G-CEM Capsule </p> <p>Self-adhesive resin cement</p> <ul style="list-style-type: none"> • Efficient self-cure mode • Easy-to-use <p>Page 10</p>	
<p>Fuji PLUS </p> <p>Resin-modified glass ionomer (RMGI) cement</p> <ul style="list-style-type: none"> • Long-term clinical success • Worry-free cement <p>Page 8</p>		<p>G-CEM LinkForce </p> <p>Adhesive resin cement</p> <ul style="list-style-type: none"> • High adhesion • High aesthetics <p>Page 11</p>	
<p>FujiCEM 2 </p> <p>Resin-modified glass ionomer (RMGI) cement</p> <ul style="list-style-type: none"> • Easy-to-use • Versatile <p>Page 9</p>		<p>G-ænial Universal Flo </p> <p>Injectable restorative composite for bonding of veneers (light-cure)</p> <ul style="list-style-type: none"> • High aesthetics • Exceptional physical properties <p>Page 12</p>	

Our solutions for temporary cementation

Zinc oxide cement	Conventional glass ionomer cement
<p>Freegenol </p> <p>Page 13</p>	<p>Fuji TEMP LT </p> <p>Page 13</p>
	

GC Fuji I

Conventional Glass Ionomer cement



Cost-efficient solution



Resin-free cement for patients allergic to methacrylate (resin)



Good solution for cementation on implant abutments. Easy removal of excesses considerably reduces the risk of peri-implantitis.

Source: Dr. Emmanuel d'Incau, France

The biocompatible, resin-free luting solution with 25 years of great clinical success.

Discover why Fuji I is still a material of choice.

Key clinical advantages

- **Chemical adhesion to tooth** structure, ensuring an excellent **marginal seal**
- **Moisture-tolerant solution** when a dry field cannot be assured
- **Quick and easy** without the need for etching and bonding steps of the tooth
- **Resin-free cement** for patients allergic to methacrylate (resin)
- **Virtually no post-operative sensitivity:** a pulp friendly solution

Practical advantages

- **Cost efficient** and **easy-to-use** solution for metal-based-restorations
- Economic consideration: €€€€€€

Clinical applications

	Metal / PFM	Zirconia / Alumina
Inlays / onlays	●	N.A.
Crowns / bridges	●	●
Posts/ inlay cores	●	-

“Fuji I is perfectly indicated for the permanent cementation of fixed elements on implant abutments. The removal of excess is simple and easy to control, thereby reducing the risk of peri-implantitis.”

DR. EMMANUEL D'INCAU, FRANCE

GC Fuji PLUS

Resin-Modified Glass Ionomer cement



Long-term clinical success



Perfect balance between protection of the tooth and adhesion



Easy removal of excess



Predictable and worry-free solution

Source: Dr. Emmanuel d'Incau, France

Efficient and proven over 20 years - a predictable, worry-free solution.

GC Fuji PLUS is one of the most well-known products and it is considered a benchmark material in its class, offering the tooth protection of a glass ionomer cement balanced with the strength and aesthetics of a resin cement.

Key clinical advantages

- **Chemical adhesion to tooth** structure, ensuring an excellent marginal seal
- **Moisture-tolerant solution** when a dry field cannot be assured
- **Quick and easy:** no need for etching and bonding of the tooth
- Rubbery consistency for **easy removal of excess**
- Virtually **no post-operative sensitivity:** a pulp friendly solution

Practical advantages

- **Proven and predictable performance** for most routine cases
- Available in different shades
- Also available in extended working time (EWT) version for long-span bridges
- Economic consideration: €€€€€€

Clinical applications

	Metal / PFM	Feldspathic / leucite ceramics*	Lithium disilicate	Zirconia / Alumina	Composite / Fibre-reinforced resin
Inlays / onlays	●	● (inlays)	●	N.A.	●
Crowns / bridges	●	-	●	●	●
Posts/ inlay-cores	●	N.A.	N.A.	●	●

"Our office has been using Fuji Plus for almost 2 decades. I don't remember any failure yet and we have cemented more than 5000 crowns & bridges!"

DR. HAN HEESSEN, NETHERLANDS

* Feldspathic & leucite-reinforced ceramics

GC FujiCEM 2

Resin-Modified Glass Ionomer cement

Simplicity is innovation! The first resin-modified glass ionomer cement in paste-paste delivery form. This easy-to-use and versatile luting cement is your perfect partner for most of your routine cementation.



Handmix

Automix



Easy dispensing

Easy removal of excess

Source: Dr. Frédéric Raux, France

Key clinical advantages

- **Chemical adhesion to tooth** structure, ensuring an excellent **marginal seal**
- **Moisture-tolerant:** a "must-have" when moisture control cannot be guaranteed
- **Quick and easy:** no need for etching and bonding of the tooth
- Rubbery consistency for **easy removal of excess**
- Virtually **no post-operative sensitivity:** a pulp friendly solution

Practical advantages

- Choose from handmix or automix: the choice is yours!
- Economic consideration: €€€€€€

Clinical applications

	Metal / PFM	Feldspathic / leucite ceramics*	Lithium disilicate	Zirconia / Alumina	Composite / Fibre-reinforced resin
Inlays / onlays	●	● (inlays)	●	N.A.	●
Crowns / bridges	●	-	●	●	●
Posts/ inlay-cores	●	N.A.	N.A.	●	●

"Why choose a very technique-sensitive bonding procedure when you can have the fantastic chemical bond to dentin, marginal seal and strength of a resin-modified glass ionomer cement? FujiCEM 2 is easy, fast and reliable."

DR. JAVIER TAPIA, SPAIN

* Feldspathic & leucite-reinforced ceramics

GC G-CEM LinkAce

Self-adhesive resin cement



Can be stored at room temperature



No need for etching and bonding of the tooth



Efficient self-cure mode for optimal polymerisation



Long-lasting aesthetic results

Source: Dr. Frédéric Raux, France

Looking for the best results even when light cannot pass through the restoration?

Choose our G-CEM LinkAce featuring highly efficient self-cure for predictable results every time! Thanks to its unique initiator and filler particles, it offers reliable bonding, long-term colour stability and high wear resistance.

Key clinical advantages

- **Efficient self-cure mode** guarantees **optimal polymerization** even under opaque or thick restorations
- **No pre-treatment** of tooth or **application of primers** for metal or zirconia
- Thin film thickness (3µm) ensures an **excellent fit of the prosthesis**
- **Excellent colour stability** for long-lasting aesthetic results

Practical advantages

- **Easy-to-use** resin cement, fulfilling most of your daily aesthetic challenges
- **Refrigeration not needed**, store at room temperature
- Economic consideration: €€€€€€

Clinical applications

	Metal / PFM	Feldspathic / leucite ceramics*	Lithium disilicate	Zirconia / Alumina	Composite / Fibre-reinforced resin
Inlays / onlays	●	●	●	N.A.	●
Crowns / bridges	●	●	●	●	●
Posts/ inlay-cores	●	N.A.	N.A.	●	●

"I use this cement most of the time for zirconia and ceramic restorations. I particularly love the fact that it is a dual-cure material with an excellent self-cure feature, ensuring high polymerisation even when transmission of light is limited. The material also offers aesthetic results and is easy to use!"

DR. ANJA BARABA, CROATIA


* Feldspathic & leucite-reinforced ceramics

GC G-CEM LinkForce

Dual-cure adhesive resin cement

Strength and aesthetics in one system for all indications, all substrates

G-CEM LinkForce is the universal adhesive resin cement that is ideal to use whenever additional retention is needed and a must-have for all CAD/CAM ceramic and hybrid ceramic blocks such as CERASMART.



Easy application of G-Premio BOND with or without prior etching

Light-curing of the thin bonding (3µm) for optimal adhesion



One universal primer for stable adhesion to all substrates

Easy seating and perfect adaptation thanks to a very low film thickness

Source: Dr. Antonio Saiz-Pardo, Spain

Key clinical advantages

- Secure adhesion in all situations with **only one system, three base elements**:
 - G-Premio BOND bonds to ALL Preparations with no compromises
 - G-Multi Primer ensures a stable adhesion to ALL Restorations
 - G-CEM LinkForce provides a strong link in ALL Indications
- Optimal light-curing feature for veneer indications
- **Efficient self-cure mode**, useful for luting more opaque or thick restorations

Practical advantages

- The universal, aesthetic and powerful solution to all your adhesive cementation challenges
- Choice of **four aesthetic shades** and their corresponding **try-in pastes**
- Economic consideration: €€€€€€

Clinical applications

	Metal / PFM	Feldspathic / leucite ceramics*	Lithium disilicate	Zirconia / Alumina	Composite / Fibre-reinforced resin	Hybrid ceramics (e.g. CERASMART)
Inlays / onlays	●	●	●	N.A.	●	●
Crowns / bridges	●	●	●	●	●	●
Veneers	N.A.	●	●	N.A.	●	●
Posts/ inlay-cores	●	N.A.	N.A.	●	●	N.A.

"I use G-CEM LinkForce as a first choice in all clinical indications. Day after day it still surprises me."

DR. CASTI, ITALY


* Feldspathic & leucite-reinforced ceramics

GC G-ænial Universal Flo

Light-cure adhesive resin cement

A versatile, injectable restorative composite for the bonding of veneers

G-ænial Universal Flo generally used as a restorative and liner material, can also be used for the bonding of veneers and some inlays and onlays where light-curing is possible.



Easy dispensing

Thixotropic consistency: easy removal of excess



Light-curing composite: high aesthetics over time

Source: Prof. Jean-François Lasserre, France

Key clinical advantages

- Injectable composite with excellent thixotropic consistency and **initial high viscosity** for **easy placement of veneers**
- Thanks to direct dispensing, the risk of bubbles at margins is limited
- High wear resistance and excellent polishability for **long-term aesthetics**
- **Good alternative to heated composite** with less risk of failure for thin indirect restorations

Practical advantages

- Streamline your inventory with this multi-purpose product!
- Broad range of shades: 15 shades

Clinical applications

	Glass-ceramics*	Hybrid ceramics (e.g. CERASMART)	Composite
Veneers	●	●	●
Inlays / onlays	● (< 2.0 mm thickness)	● (< 2.0 mm thickness)	● (< 2.0 mm thickness)

“Well-known for its optical and mechanical properties, G-ænial Universal Flo is a good solution for direct restorations and cementation (veneers, inlays and onlays). Nowadays, the practitioner is looking for a multi-purpose material and G-ænial Universal Flo fulfils this requirement. Its high viscosity allows for easy control of the excess and easy placement and positioning of the restoration.”

DR. STEPHEN KOUBI, FRANCE

* Feldspathic, leucite-reinforced ceramics and lithium disilicate

GC's temporary cements...

GC Freegenol Zinc oxide cement



Key features

- All-purpose temporary solution
- Adjustable consistency and setting times
- Easy removal of excess
- Eugenol free, avoid interaction with resin cements

GC Fuji TEMP LT

Conventional glass ionomer cement in Paste-pak cartridge



Key features

- Assured retention and easy removal
- Easy seating thanks to low film thickness (6µm)
- Rubbery consistency for easy removal of excess
- Radiopaque
- Ideal for cementation on implant abutments

Also discover GC's Primers...

Whether you prefer a specialised or a universal primer, GC has a solution.

Universal solution G-Multi PRIMER

One universal primer for **ALL** your indirect restorations.



Specialised solutions Ceramic Primer II

A specialised primer for all **aesthetic indirect restorations**. Compatible with glass-ceramic, hybrid ceramic, zirconia, alumina and composite restorations.










Metal Primer Z

The hero for all **your metal and porcelain fused to metal restorations**.













GC Luting navigation chart. Find the right cement(s) to master your luting challenges!

	Inlays & onlays				Crowns & bridges				Posts & inlay-cores			Veneers
	Metal	<ul style="list-style-type: none"> Feldspathic ceramics Leucite-reinforced ceramics 	<ul style="list-style-type: none"> Lithium disilicate Composite (e.g. GC GRADIA PLUS) 	Hybrid ceramics (e.g. CERASMART)	<ul style="list-style-type: none"> Metal (e.g. GC Initial Cast NP) Zirconia (e.g. GC Initial Zirconia disk) Alumina 	<ul style="list-style-type: none"> Lithium disilicate Composite (e.g. GC GRADIA PLUS) 	<ul style="list-style-type: none"> Feldspathic ceramics Leucite-reinforced ceramics 	Hybrid ceramics (e.g. CERASMART)	Fibre reinforced (e.g. EverStick POST & GC Fiber post)	Metal	Zirconia (e.g. GC Initial Zirconia disk)	<ul style="list-style-type: none"> Feldspathic ceramics Leucite reinforced ceramics Hybrid ceramics (e.g. CERASMART) Composite
Fuji I Conventional glass ionomer cement 	●	-	-	-	●	-	-	-	-	●	-	-
Fuji PLUS Resin-modified glass ionomer cement 	●	●	●	-	●	●	-	-	●	●	●	-
FujiCEM 2 Resin-modified glass ionomer cement 	-	(inlays)	-	-	-	-	-	-	-	-	-	-
G-CEM LinkAce Self-adhesive resin cement 	●	●	●	-	●	●	●	-	●	●	●	-
G-CEM Self-adhesive resin cement 	●	●	●	-	●	●	●	-	●	●	●	-
G-CEM LinkForce Adhesive resin cement 	●	●	●	●	●	●	●	●	●	●	●	●
G-ænial Universal Flo Adhesive resin cement 	-	●	●	●	-	-	-	-	-	-	-	●

(if less than 2.0 mm thickness)*1

Discover more features about GC luting cements

	Conventional glass ionomer		Resin-modified glass ionomer				Self-adhesive resin		Adhesive resin	
	Fuji I 	Fuji I 	Fuji PLUS 	Fuji PLUS EWT*2 	Fuji PLUS 	FujiCEM 2 	G-CEM 	G-CEM LinkAce 	G-CEM LinkForce 	G-ænial Universal Flo 
Shades	Light yellow	Light yellow	Yellow, Translucent & A3	A3	Yellow & A3	Light yellow	A2, Translucent, AO3 & BO1	A2, Translucent, AO3 & BO1	A2, Translucent, Opaque & Bleach	A1, A2, A3, A3.5, A4, BW, B1, B2, B3, C3, CV, AO2, AO3, JE & AE
Mixing time	20 sec	10 sec	20 sec	20 sec	10 sec	10-15 sec	10 sec	10-15 sec	10-15 sec	N.A.
Working time	2 min	2 min 15	2 min 30	3 min 30	2 min	2 min 15	2 min 30	2 min 50	3 min	N.A.
Final finishing time (commencing time)	4 min 30	4 min 30	4 min 30	5 min 45	4 min 15	4 min 30	4 min	4 min	4 min	N.A.

Notes: *1 For indirect restorations with thin thickness (≤ 2mm) that enable the use of a purely light-curing technique. *2: EWT: Extended WorkingTime: perfect for the cementation of bridges.

Let's Join Forces! Discover our key luting solutions for all your luting challenges

FujiCEM 2 & Fuji PLUS

Resin-modified glass ionomer cement

Proven to be a benchmark in their class. They also are ideal worry-free solutions for your routine procedures, especially for metal, porcelain fused to metal & zirconia. They are your 'must-haves', especially when moisture control becomes a challenge!



- 1 | Perfect solution for subgingival margins!
- 2 | Easy dispensing
- 3 | Easy removal of excess
- 4 | Full protection for the tooth




Source: Prof. Jean-François Lasserre (1), Dr. Frédéric Raux (2 and 4) & Dr. Emmanuel d'Incau (3), France

G-CEM LinkAce

Self-adhesive resin cement

With probably the best self-curing ability in its class, excellent polymerisation is assured, even when light transmission through the restoration is limited.



- 1 | Can be stored at room temperature
- 2 | No need for etching and bonding of the tooth
- 3 | Efficient self-cure mode for optimal polymerisation even when light transmission is limited
- 4 | Aesthetic results in time



Source: Dr. Frédéric Raux, France

G-CEM LinkForce

Adhesive resin cement





When uncompromised adhesion is needed, this product is THE solution. This universal adhesive resin cement allows you to work in a standardised way and offers a strong bond in all indications.







- 1 | Easy application of G-Premio BOND
- 2 | Light-curing of the thin bonding layer for optimal adhesion
- 3 | Universal Primer for stable adhesion to all substrates
- 4 | Easy seating and perfect adaption



Source: Dr. Antonio Saiz-Pardo, Spain

Fuji I	
	000136 Intro Package (1x Powder Light Yellow 35 g & 1x Liquid 20 ml)
	000134 Refill- 35 g Powder Light Yellow
	000135 Refill- 20 ml (25 g) Liquid
	000240 50 Capsules Light Yellow
	Fuji PLUS
	001409 Intro Package - Full Set A3 (15 g Powder A3, 7 ml Liquid & 6,5 ml conditioner)
	000231 Refill- EWT (extended working time version) 15 g Powder A3
	000219 Refill- 15 g Powder Yellow
	001410 Refill- 15 g Powder Translucent
	001411 Refill- 15 g Powder A3
	000220 Refill- 7 ml (8 g) Liquid
	000222 50 Capsules Yellow
	001418 50 Capsules A3
	
900897 Refill - 2x Paste Pak cartridge (13,3 g / 7.2 mL)	
900898 Automix Package - 2x Paste Pak cartridges (13,3g) & 44 mixing tips	
900902 20x Mixing Tips SL, Refill Pack	
900903 10x Mixing Tips SL for Endo, Refill Pack	
001573 GC Paste Pak Dispenser	
004647 GC FujiCEM 2 Dispenser	
	G-CEM
	002588 50 Capsules Assorted (20 each of A2 and Translucent, 5 each of AO3 and BO1)
	002584 50 Capsules A2
	002585 50 Capsules AO3
	002586 50 Capsules Translucent
	002587 50 Capsules BO1
004414 50x GC Capsule elongation tips	
	G-CEM LinkAce
	004857 Refill- 2 x 4.6 g (2.7 ml) A2
	004858 Refill- 2x 4.6 g (2.7 ml) Translucent
	004859 Refill- 2 x 4.6 g (2.7 ml) AO3,
	004860 Refill- 2 x 4.6 g (2.7 ml) BO1
	900668 10x GC Automix Tip Regular
	900495 10x GC Automix Tip for Endo

G-CEM LinkForce	
	009540 System Kit: 1x G-CEM LinkForce A2, 1x G-CEM LinkForce Translucent, G-Premio BOND, G-Premio BOND DCA, G-Multi PRIMER, G-CEM LinkForce Try-in Pastes (A2 & Translucent), GC Etchant & accessories
	009542 Starter Kit A2: 1x G-CEM LinkForce A2 (8.7g), G-Premio BOND, G-Multi PRIMER & accessories
	009543 Starter Kit Translucent (TR): 1x G-CEM LinkForce Translucent (8.7g), G-Premio BOND, G-Multi PRIMER & accessories
	009544 Refill A2- 1x 8.7g
	009545 Refill Translucent (TR)- 1x 8.7g
	009546 Refill Opaque- 1x 8.7g
	009547 Refill Bleach - 1x 8.7g
	009548 1x Try-In Paste 1.5 g A2
	009549 1x Try-In Paste 1.5 g Translucent (TR)
	009550 1x Try-In Paste 1.5 g Opaque
009551 1x Try-In Paste 1.5 g Bleach	
009036 G-Premio BOND 5 ml	
009552 G-Premio BOND DCA 3 ml	
	G-aenial Universal Flo
	004213 1x 2 ml (3.4 g) syringe - BW (bleach)
	004202 1x 2 ml (3.4 g) syringe - A1
	004203 1x 2 ml (3.4 g) syringe - A2
	004204 1x 2 ml (3.4 g) syringe - A3
	004212 1x 2 ml (3.4 g) syringe - AO3
	004215 1x 2 ml (3.4 g) syringe - JE (junior enamel)
	004634 30x Dispensing Tip III plastic
	004635 30x Dispensing Tip III Needle metal
	Also available in A3.5, A4, B1, B2, B3, C3, CV, AO2, AO3 & AE (adult enamel)
	Temporary Cements
	Freegenol 000087 1-1 Pack tubes (55 g Base, 20 g Accelerator, 2.5 ml Cleaner)
Fuji TEMP LT 004253 2x 13.3 g (7.2 ml) Paste Pak cartridges	
	Primers
	G-Multi PRIMER 009553 G-Multi PRIMER, 5 ml Liquid
	Ceramic Primer II 008551 GC Ceramic Primer II, 3ml Liquid
	008620 GC Ceramic Primer II Unit dose pack (10 pcs)
Metal Primer Z 009295 GC METAL PRIMER Z, 5ml Liquid	

GC EUROPE N.V. | Head Office
 Researchpark Haasrode-Leuven 1240
 Interleuvenlaan 33
 B-3001 Leuven

Tel. +32.16.74.10.00
 Fax. +32.16.40.02.14
 info@gceurope.com
 www.gceurope.com

GC UNITED KINGDOM Ltd.
 Coopers Court
 Newport Pagnell
 UK-Bucks. MK16 8JS

Tel. +44.1908.218.999
 Fax. +44.1908.218.900
 info@uk.gceurope.com
 http://uk.gceurope.com

