


















GC Luting navigation chart. Find the right cement(s) to master your luting challenges!

	Inlays & onlays				Crowns & bridges				Posts & inlay-cores			Veneers
	Metal	<ul style="list-style-type: none"> Feldspathic ceramics Leucite-reinforced ceramics 	<ul style="list-style-type: none"> Lithium disilicate Composite (e.g. GC GRADIA PLUS) 	Hybrid ceramics (e.g. CERASMART)	<ul style="list-style-type: none"> Metal (e.g. GC Initial Cast NP) Zirconia (e.g. GC Initial Zirconia disk) Alumina 	<ul style="list-style-type: none"> Lithium disilicate Composite (e.g. GC GRADIA PLUS) 	<ul style="list-style-type: none"> Feldspathic ceramics Leucite-reinforced ceramics 	Hybrid ceramics (e.g. CERASMART)	Fibre reinforced (e.g. EverStick POST & GC Fiber post)	Metal	Zirconia (e.g. GC Initial Zirconia disk)	<ul style="list-style-type: none"> Feldspathic ceramics Leucite reinforced ceramics Hybrid ceramics (e.g. CERASMART) Composite
Fuji I Conventional glass ionomer cement 	●	-	-	-	●	-	-	-	-	●	-	-
Fuji PLUS Resin-modified glass ionomer cement 	●	-	-	-	●	●	-	-	●	●	●	-
FujiCEM 2 Resin-modified glass ionomer cement 	-	(inlays)	●	-	●	●	-	-	●	●	●	-
G-CEM LinkAce Self-adhesive resin cement 	●	●	●	-	●	●	●	-	●	●	●	-
G-CEM Self-adhesive resin cement 	●	●	●	-	●	●	●	-	●	●	●	-
G-CEM LinkForce Adhesive resin cement 	●	●	●	●	●	●	●	●	●	●	●	●
G-ænial Universal Flo Adhesive resin cement 	-	●	●	●	-	-	-	-	-	-	-	●

(if less than 2.0 mm thickness)*1

Discover more features about GC luting cements

	Conventional glass ionomer		Resin-modified glass ionomer				Self-adhesive resin		Adhesive resin	
	Fuji I 	Fuji I 	Fuji PLUS 	Fuji PLUS EWT*2 	Fuji PLUS 	FujiCEM 2 	G-CEM 	G-CEM LinkAce 	G-CEM LinkForce 	G-ænial Universal Flo 
Shades	Light yellow	Light yellow	Yellow, Translucent & A3	A3	Yellow & A3	Light yellow	A2, Translucent, AO3 & BO1	A2, Translucent, AO3 & BO1	A2, Translucent, Opaque & Bleach	A1, A2, A3, A3.5, A4, BW, B1, B2, B3, C3, CV, AO2, AO3, JE & AE
Mixing time	20 sec	10 sec	20 sec	20 sec	10 sec	10-15 sec	10 sec	10-15 sec	10-15 sec	N.A.
Working time	2 min	2 min 15	2 min 30	3 min 30	2 min	2 min 15	2 min 30	2 min 50	3 min	N.A.
Final finishing time (commencing time)	4 min 30	4 min 30	4 min 30	5 min 45	4 min 15	4 min 30	4 min	4 min	4 min	N.A.

Notes: *1 For indirect restorations with thin thickness (≤ 2mm) that enable the use of a purely light-curing technique. *2: EWT: Extended WorkingTime: perfect for the cementation of bridges.