GC Fuji VIII GP

RADIOPAQUE REINFORCED GLASS IONOMER ANTERIOR RESTORATIVE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Restoration of root surface lesions.
- Restoration of Class III and V lesions.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g / g)	3,4 / 1,0
Mixing Time (sec.)	25"-30"
Working Time (min., sec.)	1' 30"
Net Setting Time (min., sec.)	2' 10"
Final Finishing Commencing Time (37°C, after starting mix) (min., sec.)	6' 00"

1. PREPARATION OF THE LESION

- Prepare the lesion using a standard Minimum Intervention approach. Mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
- Apply CAVITY CONDITIONER (GC) (10 seconds) or DENTIN CONDITIONER (GC) (20 seconds) to the bonding surfaces using a cotton pellet or sponge (Fig. 1).
- Rinse thoroughly with water. Blot away excess water with a cotton pellet or dry by gently blowing with an air syringe (Fig. 2). DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).

2. POWDER AND LIQUID DISPENSING

- The standard powder to liquid ratio is 3.4g/1.0g (1 level scoop of powder to 1 drop of liquid).
- For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
- Make sure the liquid nozzle is clean and dry before dispensing liquid. Turn the liquid bottle horizontal and hold in this position briefly to remove air bubbles. Then invert and hold the liquid bottle vertically and squeeze gently to dispense a bubble free drop of liquid.
- Close bottle caps tightly immediately after use.

3. MIXING

Dispense the required amounts of powder and liquid onto the pad. Using the plastic spatula, divide the powder into 2 equal parts. Mix the first portion with all the liquid for 10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-20 seconds (Total mix not more than 30 seconds) (Fig. 3). Working time is 1 minute 30 seconds from the start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.

4. PLACEMENT

- Transfer the mixed material to the prepared lesion using a syringe or other suitable instrument. Slightly overfill it. Avoid the incorporation of air bubbles (Fig. 4).
- Form the preliminary contour and cover with a matrix.
 - Note
 - In case no matrix is used, overfill cavity slightly.
- When set (2 minutes 10 seconds), remove the matrix and immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.) (Fig. 5).

5. FINISHING

- Final finishing under water spray using standard techniques can begin approximately 6 minutes after starting the mix (Fig. 6).
- Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.) (Fig. 7).

STORAGE

Store in a cool and dark place (4-25°C) (39,2-77,0°F). (Shelf life : powder 3 years, liquid 2 years from date of manufacture).

SHADES

A2, A3, A3.5, B2, B3, C4
Shade numbers according to Vita® shade guide.
*Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- 1-1 package : 15g powder (1), 8g (6.8mL) liquid (1), powder scoop, mixing pad (No. 22), plastic spatula.
- Bottle of 15g powder with scoop.
- Bottle of 8g (6.8mL) liquid.

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with any other glass ionomer components.
- This product should not be used with eugenol containing substances in the same cavity. It may inhibit the curing and bonding of the product.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com>
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERGE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Last revised : 10/2016

CE 0086

'GC'

MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Renaatpark Hasarode-Luven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium TEL. +32 16 74 10 00

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