Prior to use, carefully read the instructions for use.

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GC Fuji IX GP-EXTRA

RADIOPAQUE POSTERIOR GLASS IONOMER RESTORATIVE CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Class I and II restorations in deciduous teeth.
- Non-load bearing Class I and Class II restorations in permanent teeth.
- Intermediate restorative and base material for heavy stress situation in Class I and Class II cavities using sandwich laminate technique.
- Class V and root surface restorations.
- Core build-up.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	3.4 / 1.0
Mixing Time (sec.)	25-30"
Working Time (min., sec.)	2'00"
Net Setting Time (min., sec.)	2'20"
Final Finishing Commencing Time (37°C, after starting mix) (min., sec.)	6'00"

Test conditions : Temperature (23 +/-1°C)
Relative humidity (50 +/-5%)

1. POWDER AND LIQUID DISPENSING

- The standard powder to liquid ratio is 3.4g/1.0g. (1 level scoop of powder to 1 drop of liquid).
- For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
- To dispense the correct amount of liquid from the bottle, turn the bottle upside down, keep the nozzle away from the mixing pad and gently press the bottle. Discard the first 1-2 drops after opening a new bottle as the dispensed amount may not be accurate.
- Close bottles tightly immediately after use.

2. MIXING

Dispense powder and liquid onto the pad. Using the plastic spatula, divide the powder into 2 equal parts. Mix the first portion with all the liquid for 10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-20 seconds (Fig. 1).

3. RESTORATIVE TECHNIQUE

- Prepare the tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
- Apply GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) to the bonding surfaces using a cotton pellet or sponge (Fig. 2). Refer to the respective instructions for use.
- Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe. DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening) (Fig. 3).
- Mix the required amount of GC Fuji IX GP-EXTRA. Working time is 2 minutes from the start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
- Transfer cement to the preparation using a syringe or other suitable instrument. Avoid air bubbles.
- Form the preliminary contour and cover with a matrix if required (Fig. 4).
- When set, immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.).

4. FINISHING

- Final finishing under water spray using standard techniques can begin at 6 minutes after starting the mix.
- Apply a final coat of G-COAT PLUS and light cure for 20 sec. (Fig. 5, 6). Alternatively, apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) Refer to the respective instructions for use.
- Instruct the patient not to apply pressure for 1 hour.

STORAGE

Recommended for optimal performance, store at room temperature, away from heat, moisture and direct sunlight (4-25°C) (39.2-77.0°F). (Shelf life : 3 years from date of manufacture).

SHADES

A1, A2, A3, A3.5, B1, B2, B3, C4
Shade numbers according to Vita® shade guide.
©Vita is a trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- 1-1 package : 15g powder, 8g (6.4mL) liquid, powder scoop, mixing pad (No. 22), plastic spatula.
- Bottle of 15g powder with scoop.
- Bottle of 8g (6.4mL) liquid.

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. Wear gloves during use.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with any other glass ionomer components.
- Do not reuse the applicators for Conditioner to other patients.

Last revised : 07/2014

MANUFACTURED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU : GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium TEL: +32 16 74 10 00

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD.
11 Tampines Concourse, #03-05 Singapore 528729
TEL: +65 6546 7588

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