initial[™]
Zirconia Disk UHT

Instruction For Use

GC

Deep chamfer

Fig. 1

Rounded Shoulder

Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6


Fig. 7

E	4-6	NO	31-33	RU	58-60
DE	7-9	PT	34-36	SK	61-63
FR	10-12	EL	37-39	SL	64-66
IT	13-15	BG	40-42	SR	67-69
NL	16-18	CS	43-45	UK	70-72
ES	19-21	HR	46-48	TR	73-75
SV	22-24	HU	49-51	LV	76-78
DA	25-27	PL	52-54	LT	79-81
FI	28-30	RO	55-57	ET	82-84

GC Initial™ Zirconia Disk UHT (Ultra High Translucent) Zirconia CAD/CAM Disk

Prior to use carefully read the Instructions for Use

INTENDED USE

GC Initial Zirconia Disk UHT is intended for use by professional dental technicians in the fabrication of custom, all-ceramic restorations. GC Initial Zirconia Disk UHT is recommended for the manufacture of single-tooth and bridgework restorations with one pontic. GC Initial Zirconia Disk UHT can be used for anterior as well as posterior restorations.

INDICATIONS FOR USE

- Anatomically reduced crown frameworks and three-unit bridge frames with one pontic in the anterior and posterior tooth region.
- Monolithic (fully anatomical) crowns and up to three-unit bridges with one pontic in the anterior and posterior tooth region.

CONTRAINDICATIONS

Similar to all full ceramic restorations, zirconium dioxide frameworks are not suitable for patients suffering from bruxism. Do not use for frameworks with more than one pontic.

TECHNICAL DATA

	Value
3 point flexural strength (MPa)	600
Thermal Expansion Coefficient 10-6xK-1 (25-500°C)	9.8
Vickers Hardness (HV)	1250

Zirconia Disk made of Zirconium Oxide
in accordance with DIN EN ISO6872:2015:
Class II, Type 1 to 4

PREPARATION DESIGN

A. Chairside

When designing restorations, the following preparation guidelines and minimum dimensions should be maintained:

All-ceramic restorations (Fig. 1 & Fig. 2)

- 1) Prepare margins with deep chamfer or rounded shoulder.
- 2) Avoid having margins in direct occlusal contact with the opposing tooth.
- 3) Antagonist contacts must be taken into consideration.

Crowns in the anterior region (Fig. 3)

- 1) Width of the shoulder/chamfer should be at least 1.0mm.
- 2) Reduce the incisal crown third – in the incisal area – by approx. 1.5mm.
- 3) Reduce the vestibular and/or oral area by approx. 1.2mm.

Crowns in the posterior region (Fig. 4)

- 1) Reduce the anatomical shape and observe the stipulated minimum thickness.
- 2) Width of the shoulder/chamfer should be at least 1.0mm.
- 3) Reduce the incisal crown third – in the occlusal area by approx. 1.5mm.

B. Labside

When designing bridges, the following minimum dimensions should be adhered to:

- Cross-section of connectors for anteriors: 12.7 mm²
- Minimum section of connectors for laterals/posteriors (for three-element bridges): 17 mm²

Framework Design

Framework design for monolithic (fully anatomical) crowns - Staining Technique

Indication		Incisal / Occlusal	Wall thickness	Pontic Width
Anterior crown	Fig. 3	1.5 – 1.2	1.0	-
Poster crown	Fig. 4	1.5	1.0	-
Anterior bridge	Fig. 5	1.5 -1.2	1.1	11
Posterior bridge	Fig. 6	1.5	1.0	9

Framework design for anatomically reduced crown frameworks - Layering Technique

In the layering technique, the following guideline regarding minimum Zirconia UHT framework thickness and maximum layer thickness have to be observed (Fig. 7)

	A	B	C	D	E	F
Overall thickness of the restoration	1.5mm	1.5mm	1.0mm	1.5mm	1.2mm	1.0mm
Minimum framework thickness	1.0mm	0.8mm	0.7mm	0.7mm	0.6mm	0.5mm
Maximum layer thickness	0.5mm	0.7mm	0.3mm	0.8mm	0.6mm	0.5mm

PRECAUTION

1. Disk is pre-sintered. Handle with care.
2. Processing: only by skilled dental technicians. Use only suitable parameters for processing ceramics.

MILLING

Please refer to the instructions of the manufacturer of the milling unit.

Use the milling program for zirconia.

After milling, check the framework for:

- discolorations
- outbreaks
- cracks

If the framework is deficient, discontinue manufacture of the restoration

COLORING

It is advised to use water- or alcohol-based coloring liquids appropriate for green state (semi-sintered) ultra high translucent zirconia.

In order to avoid too intense coloring, it is suggested to use the brush technique rather than the dipping technique.

Assure the framework is dried before the sintering process.

SINTERING

Use a suitable high-temperature sintering furnace.

	Heat	Heat	Hold	Cool
Temperature (°C/°F)	1000°C / 1832°F	1450°C / 2642°F	1450°C / 2642°F	1000°C / 1832°F
Time (hour)	2	4.5	2	1

Note: The parameters for heat-up and cool-down depend on object size as well as the sintering furnace control system.

Please refer to the operating instructions of the respective high-temperature sintering furnace.

FINISHING AND POLISHING

The hard-sintered GC Initial Zirconia Disk UHT should be grinded/contoured further only if absolutely necessary. Only water-cooled diamond tools should be used for finishing and polishing. Avoid overheating as this might cause fracture to the material. Do not apply high pressure. Polishing of the crown surface is highly recommended to avoid damage to the opposing teeth.

VENEERING AND STAINING

For higher esthetic results we recommend veneering the hard-sintered Zirconia substructures with GC Initial Zr-FS. For full contour restorations GC Initial IQ Lustre Pastes NF / GC Initial Spectrum Stains can be used.

Please follow the instructions for use of the ceramic manufacturer.

CEMENTATION

1. Preparations

- Sandblast the inner surface of the restoration ($Al_2O_3 \leq 50\mu m$, 2bar) Wash with water spray or an ultrasonic cleaner and dry.
- Condition the sandblasted surfaces with an MDP containing primer such as CERAMIC PRIMER II or G-Multi PRIMER and dry with air syringe.

2. Cementing

- Cement with an adhesive resin cement such as G-CEM LinkForce or self-adhesive resin cement such as G-CEM LinkAce.

NOTE:

Prior to using CERAMIC PRIMER II , G-Multi PRIMER, G-CEM LinkForce or G-CEM LinkAce, refer to the respective instructions for use. In case the preparation is non-retentive, an adhesive resin cement (such as G-CEM LinkForce) is preferred.

STORAGE

Recommended for optimal performance, store in a dry ambient temperature.

PACKAGE

1 piece per box, diameter 98.5mm

UHT: Ultra High Translucency in 4 different thicknesses

- GC Initial Zirconia Disk UHT 12mm/14mm/18mm/25mm - Cold Isostatic Pressed (CIP)

CAUTION

1. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
2. Processing of GC Initial Zirconia Disk UHT produces dust which can lead to eye and skin irritation, and may cause lung damage. Always ensure that the extraction unit on your milling machine is in good working order. Wear protective goggles, gloves and face mask. Personal Protective Equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
3. Do not use autoclave to sterilize the framework.

Some products referenced in the present IFU may be classified as hazardous according to GHS.

Always familiarize yourself with the safety data sheets available at:

<http://www.gceurope.com>

or for the Americas

<http://www.gcamerica.com>

They can also be obtained from your supplier.

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3737 W. 127th Street, Alsip, IL 60803 U.S.A.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
GC ASIA DENTAL PTE. LTD.
GC AUSTRALASIA DENTAL PTY. LTD.