



Aadvia Lab Scan Accessories

'GC.'

Aadva Lab Scan Accessories

Model holder low

Used for scanning plaster models, with screw fixation



Aadva Lab Scan Accessories

Model holder high

Used for scanning plaster models, with screw fixation



Aadva Lab Scan Accessories

Reference points

Need to be attached to
Model holder Low and High
and are used by the software
for aligning



Aadva Lab Scan Accessories

Single Die Holder

Used for scanning single dies
without adjacent teeth



Aadva Lab Scan Accessories

Multi Die Holder (5 units)

Used for scanning up to 5 single dies without adjacent teeth



Aadva Lab Scan Accessories

Multi Die Holder (10 units)

Used for scanning up to 10 single dies without adjacent teeth



Aadva Lab Scan Accessories

Wax-up frame

Used for copying physical designed suprastructures or existing dentures



Aadva Lab Scan Accessories

Model plate with putty

Used for scanning plaster models,
fixated with putty



Aadva Lab Scan Accessories

Antagonist holder

Clamping device for scanning
two models in occlusion



Aadva Lab Scan Accessories

Impression scan holder

Holder for scanning physical impressions



Aadva Lab Scan Accessories

Calibration set

Used for calibration of the
ALS system



Aadva Lab Scan Accessories

Scanflag

Used for automatised
implant position scanning



Aadva Lab Scan Accessories

Scanbody

Used for implant position
scanning



Aadva Lab Scan Accessories

GC case

