



The Art of Creating  
the Foundation of  
Perfect Aesthetics

# G-ænial® Flo X

from GC

Featuring high  
radiopacity even at  
low thickness,  
excellent handling properties  
& optimal durability...

... the perfect partner  
for all your posterior  
restorations!



**GC**

# G-ænial® Flo X from GC

The perfect partner for all your posterior restorations

GC presents G-ænial Flo X – an integral part of the G-ænial range of restorative materials. Because we believe that your restorations deserve the best foundation.

## Optimal viscosity for lining

G-ænial Flo X **flows easily** in the cavity and ensures an **excellent adaptation** to the cavity walls



\*\*Tetric EvoFlow is not a trademark of GC.

## High radiopacity

G-ænial Flo X allows an easy diagnosis **even at a thickness of only 0.5mm!**



Radiopacity of enamel: 190-210% AI

G-ænial Flo X restoration: 303% AI at 1mm and 284% AI at 0,5mm thickness\*

Dr A. Coryn, Belgium

## 8 shades including 2 opaques


Developed to **mask cavity bottom discolorations**, these opaques (AO2 and AO3) are the perfect foundation for all your aesthetic demands



Note the translucency difference at 1mm thickness between Standard Shade A2 and Inside Shade AO2

## Comfortable application

This unique **ergonomic syringe design** with long & thin nozzle makes it possible to reach **even the deepest of cavities**



G-ænial Flo X also exhibits the **best balance between strength\* and viscosity** in its category; this makes it the material of choice for **all traditional flowable indications!**

\* Source: GCC R&D - Data on file



Available in **8 shades** perfectly matching those of G-ænial Posterior, and **two types of tips** to suit your preferences!



G-ænial Flo X  
Syringe 3,6g (2mL)  
Dispensing Tips III Plastic (20 pcs)

	A1	900711	A4	900715	
	A2	900712	CV	900716	
	A3	900713	AO2	900717	
	A3.5	900714	AO3	900718	
	GC Dispensing Tips III Plastic (30 pcs)	004634	GC Dispensing Tips III Needle (30 pcs)	004635	

GC EUROPE N.V.  
Head Office  
Researchpark Haasrode-Leuven 1240  
Interleuvenlaan 33, B-3001 Leuven  
Tel. +32.16.74.10.00  
Fax. +32.16.40.48.32  
info@gceurope.com  
http://www.gceurope.com

